



The Israeli tourniquet:

LST - LIFE SAVING TOURNIQUET

THE ULTIMATE SOLUTION
FOR BLEEDING CONTROL

The only tourniquet capable of stopping bleeding both from limbs and junctional areas, (shoulders and upper thighs), in one simple, innovative device.



- Dual Functionality
- Proven in Clinical and Field Testing
- Compact and Lightweight
- Patented
- FDA & CE Approved



**Shoulder Bleeding
Control**



**Limb Bleeding
Control**



**Upper thigh bleeding
control**

LST - THE ULTIMATE TOOL FOR SAVING LIVES.
ANYTIME, ANYWHERE.

LST - LIFE SAVING TOURNIQUET

UNIQUE FEATURES



Dedicated Bridge: Designed to hold the Pressure Pin securely, ensuring continuous pressure on the main artery leading to the limb. Its dynamic structure keeps the pin stable during evacuation, even when the casualty is in motion.

Cut Here: Ensures smooth evacuation by eliminating excess strap.

Pressure Pin: Specifically designed to control bleeding in cases of upper limb injuries, including complete amputations, by applying pressure on the main artery leading to the limb (subclavian or femoral).

Durable Handle and Inner Strap: Tightens the tourniquet to effectively compress the artery and control bleeding.

Durable Buckle: Ensures effective pressure application.



MEDSPECS SOLUTIONS INC.

Timely solution in changing the future

Contact Us

(02) 8525-0904
(+63) 917-5907-824
(+63) 998-8552-774

medspecs.solutions@gmail.com

Visit Us

Admin II The Colonnade Residences No. 132
Legaspi Street, Legaspi Village Makati City
1229

www.medspecs-solutions.com